

Title	Measurement of two-photon excitation spectra of fluorescent proteins with nonlinear Fourier-transform spectroscopy
Authors	Hiroshi Hashimoto <sup>1,2</sup> , Keisuke Isobe <sup>1</sup> , Akira Suda <sup>1,*</sup> , Fumihiko Kannari <sup>2</sup> , Hiroyuki Kawano <sup>3</sup> , Hideaki Mizuno <sup>3</sup> , Atsushi Miyawaki <sup>3</sup> , and Katsumi Midorikawa <sup>1</sup>
Publication	10 June 2010 / Vol. 49, No. 17 / APPLIED OPTICS
Abstract	We present measurements of two-photon excitation (TPE) spectra of various fluorescent proteins with nonlinear Fourier-transform spectroscopy. By using an ultrabroadband laser pulse with a spectrum ranging from 700 to 1100 nm, the absolute TPE spectra of six typical fluorescent proteins (SeBFP, Sapphire, eGFP, eCFP, Venus, DsRed) were measured with high spectral resolution.
Laser Quantum Product	<b>venteon ultra</b>
Institute	<sup>1</sup> RIKEN Advanced Science Institute, 2-1 Hirosawa, Wako, Saitama 351-0198, Japan <sup>2</sup> Department of Electronics and Electrical Engineering, Keio University, 3-14-1 Hiyoshi, Kohoku-ku, Yokohama 223-8522, Japan <sup>3</sup> RIKEN Brain Science Institute, 2-1 Hirosawa, Wako, Saitama 351-0198, Japan