

<b>Title</b>	<b>Few-Cycle High-Repetition Rate OPCPA For Multiphoton PEEM Towards Atto-PEEM</b>
<b>Authors</b>	<b>A. Harth<sup>1</sup>, C. Guo<sup>1</sup>, E. Mårzell<sup>1</sup>, A. Losquin<sup>1</sup>, R. Svärd<sup>1</sup>, E. Lorek<sup>1</sup>, S. Lehmann<sup>1</sup>, M. Miranda<sup>1</sup>, K. A. Dick<sup>1</sup>, P. Rudawski<sup>1</sup>, C. M. Heyl<sup>1</sup>, E. W. Larsen<sup>1</sup>, J. Ahrens<sup>2</sup>, O. Prochnow<sup>2</sup>, T. Binhammer<sup>2</sup>, U. Morgner<sup>3</sup>, J. Mauritsson<sup>1</sup>, A. L'Huillier<sup>1</sup>, A. Mikkelsen<sup>1</sup>, and C. L. Arnold<sup>1</sup></b> <sup>1</sup> Laser Quantum VENTEON, Germany <sup>2</sup> VENTEON Laser Technologies GmbH, Germany <sup>3</sup> Leibniz University of Hanover, Germany
<b>Publication</b>	<a href="https://www.osapublishing.org/viewmedia.cfm?uri=UP-2016-UTu2B.4&amp;seq=0">https://www.osapublishing.org/viewmedia.cfm?uri=UP-2016-UTu2B.4&amp;seq=0</a>
<b>Abstract</b>	We present a few-cycle high-repetition rate optical parametric amplifier for multiphoton PEEM experiments on semiconductor nanowires. This parametric amplifier is also used for the generation of high-order harmonics at 200kHz for future atto-PEEM experiments.
<b>Laser Quantum Product</b>	<b>Venteon OPCPA</b>
<b>Institute</b>	OSA Publishing